



Mälardalen University
School of Innovation Design and Engineering
Västerås, Sweden

GUI USER MANUAL

13th December 2023

Contents

1. GUI overview 1

2. Scan Button 1

2.1. Laser scanning 4

2.1.1 Cylinder Scan 4

2.1.2 Halve Sphere Scan 4

2.2. Microwave measurement 5

3. Stop Button 6

4. Save Button 7

5. Load Button 7

6. Calibration Button 7

7. Manual Control 8

7.1. robotic arm movement 8

7.2. Linear actuator movement 8

8. Log 9

9. 3D reconstruction 9


10. Network viewer 10

1. GUI overview

This GUI is not dynamic adjustable. Therefore, you need a big screen for display all functions.

Run MainProgram.py for run the GUI.

It is able to perform cylindrical scan. When you use halve sphere scan, your robot will have a collision with linear actuator. Same problem in microwave scanning when you want to scan the point on the bottom.



2. Scan Button

Perform scanning based on which tab is activating in Scanning Settings when user press the scan button. When scanning, other buttons about robot arm will be disable expect stop button. The result will be save in the Scanning result table in Result and plot in the Scanning viewer.




Figure 1: Scanning viewer

Result:

Scanning Result		Scanning Position list	
	X	Y	Z
1	56.54875366568915	7.770536572522401e-15	-19.9999999999...
2	54.2991202346041	8.046037210585721e-15	-24.9999999999...
3	51.56341642228739	8.381064302302607e-15	-29.9999999999...
4	48.14222873900293	8.800038956864273e-15	-34.9999999999...
5	44.06671554252199	9.299145375959591e-15	-39.9999999999...
6	39.08137829912023	9.909673105739796e-15	-44.9999999999...
7	32.494501466275665	1.0716332868711893e-14	-49.9999999999...
8	23.73277126099707	1.1789335353800603e-14	-54.9999999999...
9	6.414956011730197	1.3910156055124591e-14	-59.9999999999...
10	-2.5087976539589505	1.500300069143116e-14	-64.9999999999...
11	52.1406044523113	21.59734551447831	-19.9999999999...

Figure 2: Scanning result table in Result




Figure 3: Scan button process. Blue color means the functions work in other thread than the main thread.

2.1. Laser scanning

When this tab is activating, the robot will preform a scanning by auto generate scanning points from values in spin boxes.

Meaning of spin boxes:


- circle_radius = Radius of the cylinder
- azimuthPoints = Number of points in the azimuth angle
- zMin = Lowest point of the cylinder
- offset = Offset in the z-axis

2.1.1 Cylinder Scan

When this tab is activating, the robot will preform a cylindrical scanning by auto generate scanning points.

Meaning of spin boxes:

- z_stepsize = Number of mm between each z-plane
- laser_angle = The angle of the end effector to point the laser (between +90 and -90)



The screenshot shows a software interface titled "Scanning Settings:". It has two tabs: "Laser measurement" (selected) and "Microwave Measurement". In the top right corner, there is a "distance:" label with a spin box set to "2,00", which is crossed out with a large red 'X'. Below the tabs, there are four spin boxes for "circle_radius:" (120,00), "azimuthPoints:" (16), "offset:" (-20,00), and "zMin:" (-60,00). At the bottom, there are two sub-tabs: "Cylinder Scan" (selected) and "Halve Sphere Scan". Under "Cylinder Scan", there are two spin boxes: "z_stepsize:" (5) and "laser_angle:" (0,00).

Figure 4: A cylindrical scanning setting


2.1.2 Halve Sphere Scan

Warning: Don't use this function when there is a linear actuator on the bottom of the object.

When this tab is activating, the robot will preform a halve spherical scanning by auto generate scanning points.

Meaning of spin boxes:

- elevationPoints = Number of points in the elevation plane



Scanning Settings:

Laser measurement | **Microwave Measurement**

distance: 2,00

circle_radius: 120,00 | azimuthPoints: 16 | offset: -20,00 | zMin: -60,00

Cylinder Scan | **Halve Sphere Scan**

elevationPoints: 5

Figure 5: A halve spherical scanning setting


2.2. Microwave measurement

Warning: When you want to scan the point on the bottom, collision with linear actuator happen.

When this tab is activating, the robot will preform a scanning by reading the values in the Scanning Position list. If there is no item in the table, then it will start microwave measurement without moving. After measurement, robot position data will be display in Robot position and surface position labels. Also, the result of the network analyser will be display in the network viewer.

Buttons:

- Input = Read the values in the x, y, z spin box. Recalculate the coordinates to a point on the surface of the object. And print the point in the table and scanning viewer. (Warning: If there are no mesh when you pressed the button, the GUI will be crashed. Make sure you start a laser scanning or load a CSV file before the action.)
- Clear all = Remove all items in the Scanning Position list table.
- distance spin box = how far between the surface of object and the antenna in mm.



Scanning Settings:

Laser measurement | **Microwave Measurement**

distance: 2,00

X: 120,00 | Y: 120,00 | Z: -30,00

Input

Clear all

Figure 6: When microwave measurement tab is activating

- Even you perform stop action, the connection between Raspberry Pi and Yumi controller is not disconnect. Therefore, you need to restart the GUI program to disconnect. Therefore, you are not able to use any function that move the Yumi after you pressed stop button.




Figure 9: Stop button process. Blue color means the functions work in other thread than the main thread

4. Save Button

Save items in the Scanning Result table in Result as a csv file.




Figure 10: Save button process

5. Load Button

Read data from a csv file and insert to the Scanning Result table. It will also plot the point cloud base on the data and perform a 3d reconstruction.




Figure 11: Load button process

6. Calibration Button

Warning: It cause problem when robotic arm have too many cables.

Perform a calibration scanning and update the quaternions which are used for scan the object.

Note:

- The standard quaternion is [9.99954527e-01, 9.41712207e-03, 1.50357889e-03, 9.45543129e-06] in ClassGUI.py. This number is save in a csv file called quaternions.csv in GUI folder. After calibration, the value will be updated.




Figure 12: Calibration button process. Blue color means the functions work in other thread than the main thread.

7. Manual Control

Functions for control the robotic arm and linear actuator manually.

7.1. robotic arm movement

The principal of robotic arm movement is to able perform small movement ± 1 mm. But it always go back to the start position and then go to the specific position. Therefore, it is not good to use it. Use manual input scanning points instead. It also use same distance spin box same as the manual input scanning points.




Figure 13: Manual robot controller




Figure 14: Robot control buttons process. Blue color means the functions work in other thread than the main thread.

7.2. Linear actuator movement

Functions for control the linear actuator. You are able to move the linear actuator ± 1 mm or go to specific position by use Move button. Also, before you use the controller, you need to move the linear actuator down to bottom (the zero position), and then press calibration button in the tab to let the program know the linear actuator is on the zero position.

Buttons:

- Move = move to specific position using input in the spin box.
- Calibration = set the current value to 0
- Up = move linear actuator up 1mm
- Down = move linear actuator down 1mm

Manual Controller:

Robotic Arm movement

Linear actuator movement

Move to: ▲ ▼ Move **(mm)**

position: 0 Calibration

Move 1 mm: Up Down

Figure 15: Linear actuator




Figure 16: Linear actuator buttons process. Blue color means the functions work in other thread than the main thread.

8. Log

A text box that use for display messages.

9. 3D reconstruction

Because Raspberry pi does not support the draw function in open3D. Therefore, this functions is not available. But other functions like create mesh works. You can save the mesh to a .obj file by change the code in ClassGUI.py. And if you do, the 3D reconstruction data saves in a file called meshNew.obj.

10. Network viewer

A viewer for display the result of network analyser. It is also automatically save the result in the mw_data folder. Which result will be display base on which checkbox you had checked.

Buttons:

- plot = read a csv file in mw_data and plot it in the viewer




Figure 17: Network viewer